



First Medical Source

Reducing Cost of Healthcare Worldwide

Pain Management Catalog



To place an order, email us at Sales@firstmedicalsourc.com
For additional information, email us at Info@firstmedicalsourc.com
To report a problem or file a complaint, email us at Support@firstmedicalsourc.com

First Medical Source

Reducing Cost of Healthcare Worldwide
www.firstmedicalsourc.com

Manufactured for: First Medical Source, LLC, Laguna Niguel, CA 92677

© 2008-2010 by First Medical Source, LLC, All rights reserved.
PMK-1102H C 05/2010

Indication for use

The FMS pain management system is ideal for regional or local pain management for control of acute post operative pain. The FMS pump delivers a continuous fixed rate infusion of non-narcotic local anesthetic to the surgical site through a catheter which is inserted at the site during surgery. The FMS pump is DEHP free, Latex free, safe, easy to fill, simple to use and lightweight, so that patients can keep their mobility during their recovery.

First Medical Source is proud to offer a safe and less expensive alternative to existing pain management kits in the market. Our expertise and knowledge of the supply chain enable us to control cost, quality, and delivery to support your need. Please review the available models of the FMS pain management pumps and catheter kits described below and then contact First Medical Source to order or to obtain more information.

Order your selected FMS™ Pump and Catheter Kits

Catheter Kit Options

Catheter Kit Description	Kit Order Code
2.5" (15 holes) Fixed-Hub FMS Conduction Catheter™	
Single Site with a 4.5" Tear-Away Introducer	SC25
Dual Site with a Y-Connector & two ea 4.5" Tear-Away Introducers	YC25
5" (30 holes) Fixed-Hub FMS Conduction Catheter™	
Single Site with a 4.5" Tear-Away Introducer	SC50
Dual Site with a Y-Connector and two ea 4.5" Tear-Away Introducers	YC50
10" (60 holes) Fixed-Hub FMS Conduction Catheter™	
Single Site with a 12" Tear-Away Introducer	SC100
Dual Site with a Y-Connector with two ea 12" Tear-Away Introducers	YC100
Fixed-Hub FMS Standard (3 holes) Catheter™	
Single Site with a 4.5" Tear-Away Introducer	SC
Dual with a Y-Connector & two 4.5" Tear-Away Introducers	YC

Each catheter pack includes:

- FMS Catheter with female luer cap
- Introducer
- 60 cc filling syringe
- 5 cc priming/flushing syringe
- Transparent dressing
- Adhesive strip cards
- Bifurcated set (dual packs only)

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

Introducers and Pump Options

Product Description	Order Code
FMS Single Pack Sterile Tear-Away Introducers	
8" 17 G Introducer	INT8
12" 18 G Introducer	INT12
Infusion Pumps	
275 ml volume at 5 ml per hr continuous (2.2 days)	C0050L
275 ml volume at 4 ml per hr continuous (2.8 days)	C0040L
275 ml volume at 3 ml per hr continuous (3.8 days)	C0030L
275 ml volume at 2 ml per hr continuous (5.7 days)	C0020L
275 ml volume at 1.5 ml per hr continuous (7.5 days)	C0015L
275 ml volume at 1 ml per hr continuous (11 days)	C0010L
275 ml volume at 4 ml per hr with 2 ml bolus and 60 min lockout	T4060L
275 ml volume at 2 ml per hr with 2 ml bolus and 60 min lockout	T2060L
Other models can be supplied by custom order.	TBD