



To place an order, email us at [Sales@firstmedicalsourc.com](mailto:Sales@firstmedicalsourc.com)  
For additional information, email us at [Info@firstmedicalsourc.com](mailto:Info@firstmedicalsourc.com)  
To report a problem or file a complaint, email us at [Support@firstmedicalsourc.com](mailto:Support@firstmedicalsourc.com)

## First Medical Source

Reducing Cost of Healthcare Worldwide  
[www.firstmedicalsourc.com](http://www.firstmedicalsourc.com)

Manufactured for: First Medical Source, LLC, Laguna Niguel, CA 92677

© 2009 by First Medical Source, LLC, All rights reserved.  
PMK-1102H B 11/2009



## First Medical Source

Reducing Cost of Healthcare Worldwide

### Pain Management Catalog



#### Indication for use

The FMS™ pump is ideal for regional or local pain management for control of acute post operative pain. The FMS™ pump delivers a continuous fixed rate infusion of non-narcotic local anesthetic to the surgical site through a catheter which is inserted at the site during surgery. The FMS™ pump is DEHP free, Latex free, safe, easy to fill, simple to use and lightweight, so that patients can keep their mobility during their recovery.

First Medical Source is proud to offer a safe and less expensive alternative to existing pain management kits in the market. Our expertise and knowledge of the supply chain enable us to control cost, quality, and delivery to support your needs within 3-4 weeks of the order acknowledgement.

Please review the available models of the FMS™ pain management pumps and catheter kits described below and then contact First Medical Source to order or to obtain more information.

**Order your selected FMS™ Pump and Catheter Kits**

**Catheter Kit Options**

Catheter Kit Description	Kit Order Code
<b>2.5" Fixed-Hub FMS Conduction Catheter™</b>	
Single Site with a 4.5" Tear-Away Introducer	SC25
Dual Site with a Y-Connector & two 4.5" Tear-Away Introducers	YC25
Dual Site – with one Standard Catheter, one 2.5" FMS Conduction Catheter™ & two 4.5" Tear-Away Introducers	SY25
<b>5" Fixed-Hub FMS Conduction Catheter™</b>	
Single Site with a 4.5" Tear-Away Introducer	SC50
Dual Site with a Y-Connector and two 4.5" Tear-Away Introducers	YC50
<b>10" Fixed-Hub FMS Conduction Catheter™</b>	
Single Site without Introducer	SC100
Dual Site with a Y-Connector without Introducer	YC100
<b>Fixed-Hub FMS Standard Catheter™</b>	
Single Site with a 4.5" Tear-Away Introducer	SC
Dual with a Y-Connector & two 4.5" Tear-Away Introducers	YC

Each catheter pack includes:

- FMS Catheter with luer cap
- Introducer
- 60 cc filling syringe
- 5 cc priming/flushing syringe
- Securement dressing
- Adhesive strip cards
- Bifurcated set (dual packs only)

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

**FMS™ Pump Options**

Product Code	Flow Rate mL / HR	Nominal Fill volume mL	Duration
CT-0500-50	50	50	1 hour
CT-0500-100	50	100	2 hours
CT-1000-100	100	100	1 hour
CT-2000-100	200	100	0.5 hour
CT-1750-250	175	250	0.5 hour
CT-2500-250	250	250	1 hour
CT-1000-400	100	400	4 hours
CT-2000-400	200	400	2 hours
CT-2500-500	250	500	2 hours
CT-0020-060	2	60	30 hours
CT-0005-065	0.5	65	130 hours
CT-0005-100	0.5	100	200 hours
CT-0010-100	1	100	100 hours
CT-0020-100	2	100	50 hours
CT-0050-125	5	125	25 hours
CT-0010-270	1	270	270 hours
CT-0020-270	2	270	135 hours
CT-0040-270	4	270	67.5 hours
CT-0050-270	5	270	54 hours
CT-0100-270	10	270	27 hours
CT-0060-300	6	300	50 hours
CT-0040-400	4	400	100 hours
CT-0050-400	5	400	80 hours
CT-0100-400	10	400	40 hours